

BISCO

PRODUCT SPOTLIGHT



*Where Chemistry
and Oral Health Meet*



1-800-247-3368
www.bisco.com

Rx Only

TABLE OF CONTENTS

ADHESIVES/ETCHANTS

All-Bond Universal®	1
Universal Primer™	2
Select HV® Etch w/BAC, Uni-Etch® w/BAC, Etch-37™ w/BAC	3

FLOWABLES/BASE/LINER

TheraCal LC®	4
--------------------	---

CEMENTS

eCEMENT®	5
Duo-Link Universal™	6
BisCem®	7

POST & CORE

D.T. Light-Post® Illusion™ X-RO®	8
Core-Flo™ DC & Core-Flo™ DC Lite	9

TECHNIQUE

Porcelain Veneer Cementation	10
Post Cementation & Core Build-up	11
Intraoral Repairs-Porcelain, PFM or PFZr	12
Posterior Direct Composite Placement	13
Conventional Cementation of Metal/PMF/Zirconia Crowns (Retentive Preps)	14
Bonded Cementation of Inlays/Onlays, Bridges, Crowns (Short/Tapered Preps)	15

All-Bond Universal[®]

Light-Cured Dental Adhesive

All-Bond Universal is the culmination of over 30 years of adhesive research at BISCO. As a universal adhesive it can be used with both direct and indirect restorations and is formulated to be compatible with light-, dual- and self-cured materials. The versatility of All-Bond Universal makes it an indispensable part of any dental practice.

UNIQUE BENEFITS

- Not technique sensitive; use on wet, dry, or moist tooth structure
- Impressive bond strength to ALL substrates
- Ideal chemical balance for both total- and self-etch adhesion from one bottle
- Virtually no post-operative sensitivity
- No activator required


Dentistry courtesy of Tyler Lasseigne, DDS, CDT



1. Apply Z-Prime™ Plus



2. Apply All-Bond Universal



3. Apply Duo-Link Universal™

See pages 10 and 15 to reference technique.

The First
TRULY
UNIVERSAL
Dental
Adhesive



Clinical Significance:

Other single bottle, universal adhesives may need more than 1 bottle for indirect restorations, but with ALL-Bond Universal, NO activator is required

ORDERING INFORMATION

Bottle

All-Bond Universal Standard Kit B-72020K
All-Bond Universal (6ml) B-7202P

Unit-Dose Packages

All-Bond Universal Unit-Dose (50pk) B-73050K
All-Bond Universal Unit-Dose (100pk) B-73100K



ADHESIVES



Universal Primer™ Dual-Cured Adhesive

Universal Primer is a low film thickness, dual-cured adhesive/primer designed for the dentist who prefers not to light cure the adhesive layer under indirect restorations.

Clinical Significance:

A dual-cured adhesive designed to be used without having to cure the adhesive layer under indirect restorations

UNIQUE BENEFITS

- Can be used with total-etch, self-etch, or selective-etch techniques
- Clinical evaluation confirms the ease-of-use, with virtually no post-operative sensitivity
- Low film thickness (<5µm)
- No light-curing necessary when used with dual-cured or self-cure restorative material

ORDERING INFORMATION

ACE® Universal Primer Starter Kit	M-16300K
Universal Primer Part A & B Set (6ml ea.)	B-4110P
ACE Universal Primer 4 Cartridge Pack (2ml ea.)	B-41240P

Dentistry courtesy of Ross Nash, DDS, FAACD



1. Tooth preparation



2. Apply Universal Primer



3. Apply Duo-Link Universal™



4. Final restoration

See pages 11 and 15 to reference technique.

Select HV[®] Etch w/BAC, Uni-Etch[®] w/BAC, Etch-37[™] w/ BAC

Phosphoric Acid Etchants with Benzalkonium Chloride (BAC)

SELECT HV ETCH is a 35% high-viscosity phosphoric acid etchant with BAC. It is specifically formulated for pin-point placement to etch enamel when using a self-etch adhesive.

UNI-ETCH w/BAC and ETCH-37 w/BAC are 32% and 37% semi-gel phosphoric acid etchants with BAC. They are specifically formulated to be easily applied to larger surface areas and rinsed cleanly with no residue.



Clinical Significance:


For etching tooth structure prior to bonding

UNIQUE BENEFITS

- Contains benzalkonium chloride (BAC), an antimicrobial agent. In-vitro research shows it is effective against Streptococcus mutans¹.

NOTE: Inclusion of BAC has not been shown to correlate with a reduction in secondary decay in patients. In-vivo clinical studies to evaluate the effects of BAC on oral bacteria or caries have not been performed.

Dentistry courtesy of Dr. Michael Morgan



Select HV Etch w/BAC demonstrating pin-point placement on enamel

ORDERING INFORMATION

Select HV Etch w/BAC

Bulk Syringe Kit	E-59200K
Bulk Syringe Refill (30ml)	E-59160P
4 Syringe Package (5g ea.)	E-59110P

Uni-Etch w/BAC

Bulk Syringe Kit	E-56620K
Bulk Syringe Refill (30ml)	E-56621P
4 Syringe Package (5g ea.)	E-5502EBM
Bulk Bottle (30g)	E-5637EB

Etch-37 w/BAC

Bulk Syringe Kit	E-56740K
Bulk Syringe Refill (30ml)	E-56741P
4 Syringe Package (5g ea.)	E-5503EBM
Bulk Bottle (30g)	E-5638EB

¹ Dr. Daniel Chan, University of Texas Health Science Center at San Antonio Dental School. Residual Effect of 1 and 2% Benzalkonium Chloride Incorporated into an Etchant on the Susceptibility of Actinomyces viscosus T14V. 1993



Clinical Significance:

The proprietary formulation of TheraCal LC consists of tri-calcium silicate particles in a hydrophilic monomer that provides significant calcium release* making it a uniquely stable and durable material as a liner or base

ORDERING INFORMATION

TheraCal LC 4-Syringe Pack	H-33014P
TheraCal LC Syringe (4g)	H-3301P



TheraCal LC[®]

Resin-Modified Calcium Silicate Pulp Protectant/Liner Restorations

TheraCal LC is a light-cured resin modified calcium silicate designed to protect and seal the pulp for both direct and indirect pulp capping procedures.

UNIQUE BENEFITS

- Calcium release stimulates^{1*} hydroxyapatite and secondary dentin bridge formation^{2,3}
- Alkaline pH promotes healing^{2,4}
- Significant calcium release¹ leads to a protective seal^{5,7,8}
- Protects and insulates the pulp^{5,6}
- Moisture tolerant¹ and radiopaque – can be placed under restorative materials and cements

Dentistry courtesy of Dr. Ross Nash



1. Distal-occlusal caries present on an asymptomatic mandibular first premolar



2. Incomplete excavation of caries leaving affected dentin without exposing the pulp



3. Direct placement of TheraCal LC as an indirect pulp cap on entire pulpal floor



4. 1 mm layer of TheraCal LC after light-cure of 20 seconds

See page 13 to reference technique.

* Bisco has, on file, the calcium release data for TheraCal LC.

1. Gandolfi MG, Siboni F, Prati C. Chemical-physical properties of TheraCal, a novel light-curable MTA-like material for pulp capping. International Endodontic Journal. 2012 Jun;45(6):571-9.
 2. ADA definitions for direct and indirect pulp capping at: www.ada.org/en/publications/cdt/glossary-of-dental-clinical-and-administrative-terms
 3. Apatite-forming Ability of TheraCal Pulp-Capping Material, M.G. GANDOLFI, F. SIBONI, P. TADDEI, E. MODENA, and C. PRATI J Dent Res 90 (Spec Iss A):abstract number 2520, 2011 (www.dentalresearch.org)
 4. Okabe T, Sakamoto M, Takeuchi H, Matsushima K (2006) Effects of pH on mineralization ability of human dental pulp cells. Journal of Endodontics 32, 198-201.
 5. Sangwan P; Sangwan A; Duhan J; Rohilla A. Tertiary dentinogenesis with calcium hydroxide: a review of proposed mechanisms. Int Endod J. 2013; 46(1):3-19
 6. Selcuk SAVAS, Murat S. BOTSALI, Ebru KUCUKYILMAZ, Tugrul SARI. Evaluation of temperature changes in the pulp chamber during polymerization of light-cured pulp-capping materials by using a VALO LED light curing unit at different curing distances. Dent Mater J. 2014;33(6):764-9.
 7. Cantekin K. Bond strength of different restorative materials to light-curable mineral trioxide aggregate. J Clin Pediatr Dent. 2015 Winter;39(2):143-8.
 8. Mechanical Properties of New Dental Pulp-Capping Materials Over Time. M. NIELSEN, R. VANDERWEELE, J. CASEY, and K. VANDEWALLE, USAF, JBSA-Lackland, TX, J Dent Res 93(Spec Iss A): 495, 2014 (www.dentalresearch.org)

eCEMENT®

Adhesive Cementation System for Lithium Disilicate Restorations

The eCEMENT system is indicated for luting of lithium disilicate (e.g. IPS e.max®*) veneers, inlays, onlays, crowns, bridges and all CAD/CAM blocks and pressable ingots.

UNIQUE BENEFITS

- > Universal and versatile
- > Easy removal of excess cement
- > Low film thickness
- > Radiopaque

Your One Kit For Lithium Disilicate Restorations (e.g. IPS e.max®*)



Clinical Significance:

The eCEMENT kit offers the clinician a convenient, all in one approach to cementing IPS e.max®* indirect restorations with eCEMENT light-cured and dual-cured cements, Porcelain Primer, All-Bond Universal®, Porcelain Etchant and Select HV® Etch

Courtesy of Adamo E. Notarantonio, DDS



1. Apply All-Bond Universal



2. Remove excess cement



3. Final restorations

See page 10 to reference technique.

ORDERING INFORMATION

System Kit	C-52400K
A1 LC Syringe (4.5g)	C-52010P
Translucent LC Syringe (4.5g)	C-52120P
Milky Bright LC Syringe (4.5g)	C-52230P
Universal DC Dual-Syringe (8g)	C-52340P



*IPS e.max is a registered trademark of IVOCLAR VIVADENT, Inc.
eCEMENT is a registered trademark of BISCO, Inc.
There is no sponsorship, affiliation or connection between BISCO, Inc. and IVOCLAR VIVADENT, Inc

CEMENTS



Highly radiopaque
with easy clean-up

Clinical Significance:

Extremely high degree of conversion in both light- and self-cured modes ensures a strong, long lasting restoration

Duo-Link Universal™

Resin Luting Cement

The Duo-Link Universal Kit is an adhesive cement system that contains all the components needed for cementation of indirect restorations*, including primers for restoration treatment.

UNIQUE BENEFITS

- Highly radiopaque with easy clean-up
- Ideal for all CAD/CAM restorations

ORDERING INFORMATION

System Kit with All-Bond Universal	A-19620K
System Kit with Universal Primer	A-19710K
System Kit with Ace Universal Primer	A-19730K
System Kit with Ace Universal Primer and Dispenser	M-16200K
Universal Dual-Syringe (8g)	A-19030P
Milky White Dual-Syringe (8g)	A-197MWP

Courtesy of Darren D. Simpson, DDS



1. Duo-Link Universal dispensed on the crown



2. Apply All-Bond Universal®, air dry and light cure



3. Restoration seated. Excess cement peeled away



4. Final restoration

See page 15 to reference technique.

*It is recommended to use BISCO's CHOICE™ 2 for veneer cementation.



Duo-Link Universal Kit with All-Bond Universal



BisCem[®]


Self-Adhesive Resin Cement

BisCem is a self-etching, self-adhesive, dual-cured resin luting cement exclusively formulated for luting crowns, bridges, inlays, onlays and posts.

UNIQUE BENEFITS

- A dual-cured, self-adhesive cement specifically formulated for excellent handling properties and easy clean-up
- Radiopaque
- Low film thickness ensures the restoration is completely seated
- Cost effective alternative designed to bond to a multitude of substrates, including metals, composites, porcelain and amalgam
- Glass Ionomer benefits with resin cement strength. Contains fluoride

Courtesy of Adamo E. Notarantonio, DDS



1. Fill restoration with BisCem



2. Seat the restoration and light cure



3. Remove excess cement



Clinical Significance:

No etching, priming or bonding saves time by reducing clinical steps

ORDERING INFORMATION

Translucent Dual-Syringe (8g)	D-45011P
Opaque Dual-Syringe (8g)	D-45012P

See page 14 to reference technique.



Clinical Significance:

BISCO's fiber posts have been placed in over 3,400 cases as part of a decade long multi-center clinical trial. Three (3) post fractures and zero (0) root fractures have been documented in the clinical trials ▲

D.T. Light-Post® Illusion™ X-RO®

Double-Tapered Radiopaque Translucent Fiber Post System

The D.T. LIGHT-POST ILLUSION X-RO* from BISCO maintains all the features and benefits of the clinically proven original D.T. LIGHT-POST with the added benefit of the patented** color-on-command technology and even greater radiopacity. When the post is placed, the intrinsic color of the post disappears, allowing it to blend into the natural dentition. If retreatment or removal is ever required, the color of the post will reappear on command by simply using cold water.

UNIQUE BENEFITS

- Translucency enhances esthetics and light transmission
- Double-tapered design closely conforms to the shape of the prepared canal to minimize dentin removal and cement layer
- Unidirectional, pretensed, quartz fibers bound in resin matrix gives it a greater flexural strength than metal posts
- Maximum strength – no notches or serrations to compromise strength

ORDERING INFORMATION


System Kit with Core-Flo™ DC Lite	G-91860K
System Kit with Core-Flo™ DC	G-91870K
BisCem® Kit	G-91840K
Intro Kit	G-91820K

D.T. Light-Post Illusion X-RO (10 pk)

#0.5 (Black)	X-81810P
#1 (Red)	X-81811P
#2 (Yellow)	X-81812P
#3 (Blue)	X-81813P

D.T. Drills (2 pk)

Pre-Shaping/#0.5 (Black)	X-80784P
#1 (Red)	X-80781P
#2 (Yellow)	X-80782P
#3 (Blue)	X-80783P



1. Create post space after gutta percha removal



2. Try-in and radiograph



3. Optional: Trim excess length



4. Clean, etch, and prime canal with Universal Primer™

See page 11 to reference technique.

* D.T. LIGHT-POST/D.T. LIGHT-POST ILLUSION X-RO is manufactured by RTD.

** Patents: EP 1 115 349, EP 0 432 001, EP 1 776 933, U.S. 5,328,372, U.S. 5,890,904 & Patents Pending.

▲ Data on file. BISCO, Inc

CORE-FLO™ DC / CORE-FLO™ DC Lite

Dual-Cured Core Build-Up Material

CORE-FLO DC and CORE FLO DC Lite are dual-cured, fluoride-containing core materials that are ideal for core build-ups, post cementation and as a dentin replacement material.

UNIQUE BENEFITS

- High compressive & flexural strength
- Stackable (Core-Flo DC) or self-leveling (Core-Flo DC Lite) for optimal handling
- Easy to prepare – cuts like dentin
- All-in-one procedure: cement post to core build-up
- Radiopaque



Clinical Significance:

Dual-cured, flowable or low viscosity, core materials save time when cementing posts and building the core

Dentistry courtesy of Dr. Jack D Griffin Jr.



1. Enamel etched with Select HV® Etch



2. Core-Flo DC Lite applied after use of Universal Primer™




3. Core-Flo DC Lite placed

See page 11 to reference technique.



ORDERING INFORMATION

Core-Flo DC Lite	
System Kit with Universal Primer	AB-17410K
System Kit with Ace Universal Primer and Dispenser	M-17420K
System Kit with Ace Universal Primer	AB-17412K
Natural/A1 Dual-Syringe (8g)	A-17801P
Opaque White Dual-Syringe (8g)	A-17803P
Core-Flo DC	
System Kit with Universal Primer	M-23412K
Natural/A1 Dual-Syringe (8g)	A-23011P
Opaque White Dual-Syringe (8g)	A-23012P
Blue Dual-Syringe (8g)	A-23013P

Porcelain Veneer Cementation*



OR



1. Prepare the internal surface of the restoration as instructed by the laboratory:


a. Porcelain/Lithium Disilicate:
If needed, etch the veneer with hydrofluoric acid (4% **Porcelain Etchant** or 9.5% **Porcelain Etchant**) according to the manufacturer's instructions.

b. Apply a thin coat of silane (**Porcelain Primer**) to the internal surface. Wait for 30 seconds, or dry with (warm) air.


2. Remove the temporary restorations and clean the preps with **Cavity Cleanser™** & pumice slurry. Rinse thoroughly.




3. Try in the restorations using the corresponding shade of water-soluble try-in paste (**CEMENT® Try-In** or **Choice 2™ Try-In Paste**). Remove the veneers and either clean with etchant (**Uni-Etch® w/BAC** or **Select HV® Etch w/BAC**) and rinse thoroughly, or ultrasonicate in water or alcohol for 2-3 minutes. Dry the restorations.



4. Isolate. Etch preps with **Uni-Etch w/BAC** or **Select HV Etch w/BAC** for 15 seconds, then rinse thoroughly.




5. Dispense **ACE® Universal Primer™** & mix (or **All-Bond Universal®**) in a mixing well. Apply 2 separate coats, agitating each coat for 10-15 seconds.



6. Gently air dry until there is no visible movement of the adhesive. Then thoroughly air dry with greater air pressure. The surface should appear shiny, otherwise apply additional coats.



7. Light cure for 10 seconds.



8. Apply a thin layer of HEMA-free resin (**Porcelain Bonding Resin**) to the internal surface of the veneer. Do NOT light cure.



9. Using a light-cured resin cement (**CEMENT LC** or **Choice 2**), fill the internal surface of the restoration with the cement.




10. Seat with gentle, passive pressure and light-cure for 2-3 seconds to tack the veneer into place. Remove excess cement, then light cure each veneer for 40 seconds.


* Refer to instructions for complete details.

See page 5 for more product information.

Post Cementation & Core Build-up*




1. If using the TOTAL-ETCHING TECHNIQUE, etch the canal with Uni-Etch® w/BAC or Select HV® Etch w/BAC for 15 seconds. If using the SELF-ETCHING TECHNIQUE, proceed to the next step.




2. Clean and irrigate. Rinse canal with water. Remove excess water with a brief burst of air and paper point(s) to remove any pooled water from the canal.




3. Dispense Universal Primer™. Mix adhesive.




4. With an endodontic brush, apply and scrub 2 separate coats of Universal Primer into the canal. Remove excess pooling of material with paper points and light suction. Strongly air dry.




5. Determine the post length according to the manufacturer's instructions.




6. Coat the post with Universal Primer, air dry.



7. Inject Core-Flo™ DC or Core-Flo™ DC Lite into the canal.



8. Coat the apical end of the post with Core-Flo™ DC or Core-Flo™ DC Lite. Seat the post into the canal gently.



9. Continue to express Core-Flo DC or Core-Flo DC Lite around the post and build up the core to the desired height and shape. Remove the excess material and light cure for 40 seconds.


* Refer to instructions for complete details.

See pages 8 and 9 for more product information.

Intraoral Repairs

Porcelain, PFM or PFZ_r*

1. Isolate the area to be repaired. Remove the glaze and bevel (45 degrees) the porcelain around the area to be repaired. Sandblast or abrade with a coarse diamond bur. Rinse with water and air dry.



2. Place **Barrier Gel** on the gingival tissue and porcelain to protect areas where etching is not desired.

3. Apply **9.5% Porcelain Etchant** to the dry porcelain surface for 90 seconds. Suction the etchant with high volume evacuation, then rinse with water and air dry. The etched porcelain should appear dull and frosted.

4. Clean the etched porcelain with **Uni-Etch® w/ BAC** or **Select HV® w/BAC** and agitate for 20 seconds to remove any salts. Rinse and dry thoroughly.



5. Apply 1 coat of silane (**Porcelain Primer** or **Bis-Silane™**) to the etched porcelain and allow to dwell for 30 seconds. Dry with (warm) air syringe.

6. Apply 1 coat of **Z-Prime™ Plus** to the exposed metal/zirconia/alumina and dry with an air syringe for 3-5 seconds.

7. If metal masking is required, dispense one drop each of **Catalyst & Base** onto a mixing pad and mix with a brush tip. Apply a thin coat of **Opaker** only to the metal surface. Allow to self cure, or light cure for 5 seconds.



8. Apply a thin layer of **Porcelain Bonding Resin** to the repair site. Spread a hybrid composite (**Aelite™ All-Purpose Body**) evenly over the surface and light cure, finish/polish.

9. Place **Bis-Cover™ LV** to seal and glaze the composite.

* Refer to instructions for complete details.

Visit bisco.com to learn more about BISCO's Intraoral Repair Kit.


Posterior Direct Composite Placement*

Select the desired type and shade of **Ælite™ Composite** prior to isolating and preparing the tooth:

LS Posterior: for Class I, Class II and Class VI

All-Purpose Body: Universal Composite for Class I - Class VI


Aesthetic Enamel: Enamel Layer for Class I - Class VI




1. Clean the preparation with a slurry of pumice and water or **Cavity Cleanser™**. Rinse and dry the preparation.



2. If pulp is exposed; achieve hemostasis, apply 1 mm layer of **TheraCal LC®** directly on and around the exposed pulp, then light cure for 20 seconds. If prep is near the pulp, apply a 1mm layer of TheraCal LC as a liner and light cure for 20 seconds.




3. Dispense adhesive (**All-Bond Universal®** or **ACE® All-Bond SE®**) in a mixing well. Apply 2 coats, agitating each coat for 10-15 seconds and air dry. Light cure for 10 seconds.



4. Using the selected shades, place 2 mm increments of composite (**Aelite All-Purpose Body** or **Aelite LS Posterior**) into the cavity preparation. By lightly wetting your instrument with **Modeling Resin**, all of your direct composite restorations can be easily placed and manipulated (sculpted). Light cure each increment.



5. For the final increment, place 1-2mm of enamel composite (**Aelite Aesthetic Enamel**) and light cure.




6. Check the occlusion and polish, as needed. Etch the surface; apply a surface sealant (**BisCover™ LV** or **Fortify™ Plus**) and light-cure.


* Refer to instructions for complete details.

Visit bisco.com to learn more about BISCO's Aelite Composites.


Conventional Cementation of Metal/PFM/Zirconia Crowns (retentive preps)*




1. Remove the temporary restoration and clean the prep with **Cavity Cleanser™** & pumice slurry. Rinse thoroughly.



2. Apply one coat of **Z-Prime™ Plus** to the internal surface of the restoration and air dry for 3-5 seconds.



3. Try-in the restoration; clean the internal surface with phosphoric acid etchant. Rinse thoroughly.



4. Dispense **Bis-Cem®** into the internal surfaces of the restoration. Seat with gentle, passive pressure and remove excess cement while holding restoration in place.




5. Allow the cement to self-cure, or light-cure each surface of the tooth for 20 seconds.

* Refer to instructions for complete details.

See page 7 for more product information.

Bonded Cementation of Inlays/Onlays, Bridges, Crowns (short/tapered preps)*






OR



1. Prepare the internal surface of the restoration as instructed by the laboratory:

- a. Porcelain/Lithium Disilicate: Apply a thin coat of silane (**Porcelain Primer**) to the internal surface. Wait for 30 seconds, or dry with (warm) air.
- b. Metal/Zirconia/Indirect Composite: Apply one coat of **Z-Prime™ Plus** to the internal surface of the restoration and air dry for 3-5 seconds.



2. Try-in the restoration, then clean the internal surface with (**Uni-Etch® / Select HV® Etch**) etchant. Rinse thoroughly.



3. Remove the temporary restoration and clean the prep with **Cavity Cleanser™** & pumice slurry. Rinse thoroughly.

4. If prep includes enamel, selectively etch enamel with **Select HV® Etch w/BAC** for 15 seconds, suction then rinse thoroughly.

5. Dispense **Universal Primer™** & mix (or **All-Bond Universal®**) in a mixing well. Apply 2 separate coats, agitating each coat for 10-15 seconds.



* Refer to instructions for complete details.

6. Gently air dry until there is no visible movement of the adhesive. Then thoroughly air dry with greater air pressure. The surface should appear shiny, otherwise repeat step 5. Light cure for 10 seconds.

7. Using a dual-cured resin cement (**Duo-Link Universal™** or **zCement® Dual-Cure**), fill the internal surface of the restoration and/or the prep with the cement. Seat with gentle, passive pressure and remove excess cement while holding restoration in place. Allow the cement to self-cure, or light-cure each surface of the tooth for 40 seconds.

See pages 2, 3, 5, and 6 for more product information.

Visit www.bisco.com to learn more about BISCO products.

ADHESIVES

[Universal Adhesive \(Self-Etch or Total-Etch\)](#)

All-Bond Universal[®]

Universal Primer[™]

[Self-Etch Adhesive](#)

ACE[®] All-Bond SE[®] & All-Bond SE[®]

[Total-Etch Adhesives](#)

ACE[®] All-Bond TE[™] (All-Bond 3[®])

All-Bond 2[®]

One-Step[®] & One-Step[®] Plus

[Supportive Products](#)

Aqua-Prep[™] F

Cavity Cleanser[™]

SEMI-GEL ETCHANTS

Select HV[®] Etch

Etch-37[™] w/BAC

Uni-Etch[®] w/BAC

DENTIN DESENSITIZER

BisBlock[™]

COMPOSITES

AELITE[™] Aesthetic Enamel

AELITE[™] All-Purpose Body

AELITE[™] LS Posterior

AELITEFLO[™] & AELITEFLO[™] LV

BisFil[™] 2B

Characterization Tints

Reveal[™]

[Supportive Products](#)

Fortify[™] & Fortify[™] Plus

Modeling Resin

ESTHETIC SEALANT & LIQUID POLISH

BisCover[™] LV

Where Chemistry and Oral Health Meet

INSTRUMENTS

NB Composite Instruments

FINISHING PRODUCTS

R.A.P.T.O.R.[®]
Finishing Discs

CEMENTATION

BisCem[®]
C&B[™] Cement
Choice[™] 2
Duo-Link Universal[™]
eCEMENT[®]
Resinomer[™]
[Supportive Products](#)
Z-Prime[™] Plus
Intraoral Repair Kit
Porcelain Etchant
Porcelain Primer/Bis-Silane[™]

POST & CORE

Bis-Core[™]
Core-Flo[™] DC
Core-Flo[™] DC Lite
D.T. Light-Post[®]
D.T. Light-Post[®] Illusion[™] X-RO[®]
Light-Core[™]

PROVISIONAL MATERIALS

Pro-V[®]

PULP CAPPING/LINING MATERIAL

TheraCal LC[®]

Federal law restricts these devices to sale by or on the order of a licensed Dentist.
Not all products are available in all countries.

www.bisco.com



BISCO DENTAL PRODUCTS

1100 W. Irving Park Rd.
Schaumburg, IL 60193, USA
1-800-247-3368 • www.bisco.com



Visit us at www.bisco.com for a complete product catalog.

Rx Only

MC-50161CA