CALIBRATION CARDS CHART

Calibration card «A», High value & Low value	iPhone 5S
Calibration card «B», High value & Low value	iPhone 4S, iPhone 5, Samsung Galaxy S5
Calibration card «C», High value & Low value	iPhone 4, Samsung Note 3, Samsung Note 4
Calibration card «D», High value & Low value	iPhone 6, iPhone 5S, Samsung Galaxy S4
Calibration card «E»	iPhone 6S, Samsung Galaxy S6 – S6edge

EN

You'll find in your set Smile Capture an assortment of 8 calibration cards (4x2) called «A», «B», «C» and «D» and one calibration card «E». Because each smartphone is using a different camera with a different interpretation of colors, it made sense to create different greys. You'll find here above a **chart indicating the matching between the type of card and the model of smartphone**. Please understand this is a basic information, non-exhaustive. The color of the ambient light will have an impact on the result of your picture. For example, if you feel the result of your pictures on the display is too yellow, then we suggest you to make a calibration test with a card containing more yellow in the grey – the effect will be a reduction of the yellow in your picture. **In any event, our advice is to make a few tests in order to understand the impact the grey of the calibration card will have on your pictures – this is to say we encourage you not to limit yourself to the basic recommendations.**

High Value – Low Value

For the A,B,C and D calibration cards, you are receiving one model called « high value » and another model called « low value ». In the case the ambient light in the room is rather low, we suggest to use a « low value » card. On the contrary if the ambient light is bright, you'll better use a « high value » card.

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