**EN**
**UltraSeal XT® plus™**
**Pit & Fissure Sealant and Flowable Composite and PrimalDry Drying and Priming agent**

Description:

UltraSeal XT plus is a 58% filled, light-cure, radiopaque, methacrylate-based, thixotropic resin sealant. It is used with PrimalDry drying and priming agent to aid in avoiding moisture contamination which can cause microleakage and poor retention.

Indications for Use:

Use UltraSeal XT® plus™ for prophylactic sealing of pits and fissures. It can also be used for microrestorative of “initial layer” of composite restorations.

Directions for Use:

Pit and fissure sealant

1. Apply UltraSeal™ Brush tip or Blue Micro® tip to Ultra-Etch® 35%.
2. Attach Blue Micro tip or Black Micro® FX™ to PrimalDry® syringe.

3. Apply Inspiral Brush tip to UltraSeal XT plus syringe. Adjust Inspiral Brush tip bristles to desired length by grasping and moving inward or outward.

4. Using oil-free pump, air polishing prophylaxis system, or thin diamond fissure/prep bur, clean fissures.
5. Isolate teeth being treated to avoid saliva contamination.
6. Apply teeth to fissures for 20 seconds.

7. Thoroughly rinse and dry.

NOTE: If using sodium bicarbonate in air polishing system, repeat with a 5 second etch and thoroughly rinse.

8. Apply PrimalDry for 5 seconds, clear air syringe then air dry occlusal surface. DO NOT RINSE.
NOTE: If using a bond agent, do not use PrimalDry.
9. Express a small drop of UltraSeal XT plus to the brush tip prior to applying intraorally.

10. Use a scrubbing motion to apply resin to deep fissures.

11. Cure using curing light.

12. Light cure with VALO® LED Curing Light on Standard Power for 20 seconds or 3 seconds on Xtra Power mode. (For lights with output <600mW/cm² cure 20 seconds or for lights with output <600mW/cm² cure for 40 seconds.)

12. Adjust occlusion if necessary.
Super Adhesive Initial Layer for Restorative

1. Apply dentin bonding agent according to manufacturer instructions.

NOTE: Do not use PrimalDry with a dentin bonding agent.

2. Attach Inspiral Brush tip to UltraSeal XT plus syringe.

3. Verify flow prior to applying intraorally.

4. Apply first layer of UltraSeal XT plus. 5 – 1.5 mm thick and light cure 20 sec.

5. Light cure with 2-4mm increments as needed.

6. Follow cure with VALO® LED Curing Light on Standard Power for 20 seconds or 3 seconds on Xtra Power mode. (For lights with output <600mW/cm² cure 40 seconds or for lights with output <600mW/cm² cure for 20 seconds.)

7. For detailed anatomy, apply paste composite for final layer. Carefully shape restorations with finishing burrs and polish.
Precautions and Warnings:

1. Carefully read and understand all instructions before using UltraSeal XT® plus™.

2. UltraSeal XT plus™ is not for use only.

3. Clinician and patient should wear UV protective eyewear when curing sealant.

4. Verify flow of all syringes prior to applying intraorally. If resistance is met, replace tip and re-check. Use only recommended tips.

5. For dry use, store at room temperature. Refrigerate bulk storage.

6. Refriger refrigerated products to room temperature before use.

7. Redirect overhead light to prevent premature polymerization.

8. Resins can be sensitizing. Avoid repeated contact of uncured dental resin with skin. Do not use on patients with a known sensitivity to acrylates or other resins. If allergic reaction, dermatitis or rash develops, consult a physician.

9. One application may set in tip after a few hours. To avoid clogging, remove tip and re-pace with Luer Lock caps between each use.

10. PrimalDry is flammable. Handle with care.

11. Do not use after expiration date noted on containers.

12. Re-pace, disinfect and wipe syringe with an intermediate level disinfectant between uses. If disposable syringe cure used, remove tip, re-pace, and discard syringe cure and tip.

13. Tips are disposable. To avoid cross-contamination, do not reuse tips.

14. Dispose of tips and empty syringes properly.

15. Keep products out of direct heat/sunlight.

16. Keep out of reach of children.

NOTE: For SDS and additional information about using our products, please go to [www.ultradent.com](http://www.ultradent.com)

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