## LIGHT-CURED: THIN RESTORATION WITH HIGH TRANSLUCENCY\* (<2mm Thick)



## RESTORATION PREPARATION



restoration, remove cap from PORCELAIN to the internal surface of the restoration and allow surface of the restoration (phosphoric acid or 2 ETCHANT, securely attach black plastic dispensing to dwell for 30 seconds. Dry with air syringe for minutes ultrasonic cleaning with ethanol will tip and verify etchant flow prior to application. 3-5 seconds. Apply PORCELAIN ETCHANT (4% HF) to the dry bonding surface; 25 seconds for Lithium Disilicate (or contact laboratory for recommended etch times for other ceramics). Rinse and dry.





STEP 1A: If the lab has not already etched the STEP 1B: Brush on 1 coat of PORCELAIN PRIMER STEP 1C: After try-in, decontaminate the internal remove organic contaminants while leaving the PORCELAIN PRIMER intact). Rinse and air dry.


## TOOTH PREPARATION



STEP 2A SELF-ETCH: Prepare cavity. Wash thoroughly with water spray. Use an absorbent pellet or high volume evacuation (HVE) for 1-2 seconds, to remove excess water (do not desiccate). NOTE: When bonding predominately enamel, the total-etch technique is recommended.



STEP 2A TOTAL-ETCH: Prepare cavity. Wash thoroughly with water spray. Attach blue dispensing tip. Etch enamel and dentin using SELECT HV ETCH for 15 seconds. Rinse thoroughly. Remove excess water by blotting the surface with an absorbent pellet or high volume evacuation (HVE) for 1-2 seconds, leaving the preparation visibly moist.



UNIVERSAL into a clean well. Apply two separate air-drying with an air syringe for at least 10 coats of ALL-BOND UNIVERSAL, scrubbing the seconds; there should be no visible movement of preparation with a microbrush for 10-15 seconds the adhesive. The surface should have a uniform per coat. Do not light cure between coats.



STEP 2B: Dispense 1-2 drops of ALL-BOND STEP 2C: Evaporate excess solvent by thoroughly STEP 2D: Light cure for 10 seconds. glossy appearance; otherwise, apply an additional coat of ALL-BOND UNIVERSAL and air dry.



# MENTATION



shade of **CEMENT Light-Cured Resin Cement** to seconds to tack into place. Remove the excess lingual and 40 seconds from the facial, at a distance the internal surface of the restoration. Carefully cement. place restoration and gently guide into place with passive pressure.



STEP 3A: Apply a generous amount of the selected STEP 3B: Light cure each of the restorations for 3-5 STEP 3C: Light cure for up to 40 seconds from



of approximately 0-10 mm.

\* Refer to instructions for use for a detailed technique per indication.



# **DUAL-CURED: THICKER RESTORATIONS WITH HIGHER OPACITY\***



## RESTORATION PREPARATION



restoration, remove cap from PORCELAIN ETCHANT, to the internal surface of the restoration and allow surface of the restoration (phosphoric acid or securely attach black plastic dispensing tip and to dwell for 30 seconds. Dry with air syringe for 2 minutes ultrasonic cleaning with ethanol will verify etchant flow prior to application. Apply 3-5 seconds. PORCELAIN ETCHANT (4% HF) to the dry bonding surface; 25 seconds for Lithium Disilicate (or contact laboratory for recommended etch times for other ceramics). Rinse and dry.





STEP 1A: If the lab has not already etched the STEP 1B: Brush on 1 coat of PORCELAIN PRIMER STEP 1C: After try-in, decontaminate the internal remove organic contaminants while leaving the PORCELAIN PRIMER intact). Rinse and air dry.

## TOOTH PREPARATION



STEP 2A SELF-ETCH: Prepare cavity. Wash thoroughly with water spray. Use an absorbent pellet or high volume evacuation (HVE) for 1-2 seconds, to remove excess water (do not desiccate). NOTE: When bonding predominately enamel, the total-etch technique is recommended.





STEP 2A TOTAL-ETCH: Prepare cavity. Wash thoroughly with water spray. Attach blue dispensing tip. Etch enamel and dentin using SELECT HV ETCH for 15 seconds. Rinse thoroughly. Remove excess water by blotting the surface with an absorbent pellet or high volume evacuation (HVE) for 1-2 seconds, leaving the preparation visibly moist



UNIVERSAL into a clean well. Apply two separate air-drying with an air syringe for at least 10 seconds; coats of ALL-BOND UNIVERSAL, scrubbing the there should be no visible movement of the adhesive. preparation with a microbrush for 10-15 seconds The surface should have a uniform glossy appearance; per coat. Do not light cure between coats.



STEP 2B: Dispense 1-2 drops of ALL-BOND STEP 2C: Evaporate excess solvent by thoroughly STEP 2D: Light cure for 10 seconds. otherwise, apply an additional coat of ALL-BOND UNIVERSAL and air dry.



## /ENTATION



cement directly to the tooth preparation.



STEP 3A: Apply @CEMENT Dual-Cured Resin STEP 3B: Fully seat the restoration and gently STEP 3C: Optional: after excess cement has been restoration. For inlays, it may be easier to apply the spot curing the margins for 2-3 seconds per quarter cured for up to 40 seconds. surface (mesio-facial, disto-facial, disto-liqual, mesio-lingual).



Cement to the bondable surfaces of the remove the excess cement with a brush prior to removed each surface of the restoration may be



Use Only

Consult Instructions For Use

INTC-050R2 Rev. 2/17

REF C-52410K-EU: **C € 0459** REF C-52400K: Not CE Marked

\* Refer to instructions for use for a detailed technique per indication.