

# Ultrasonic Cleaning Bath Compatible Solutions List

Updated: 11/17/15

Product	Manufacturer	Test Results
Monarch Enzymatic Cleaner	Air Techniques	Compatible
OptiSonic™ Super-Concentrated General Purpose Ultrasonic Cleaning Solution	Benco	Compatible
Benco Enzyme Tablets	Benco	Compatible
Benco General Purpose Ultrasonic Cleaning Solution	Benco	Compatible
BiozymeLT Enzyme Detergent	Biotrol	Compatible
BioSonic Enzymatic Ultrasonic Cleaner	Coltene/Whaledent Inc	Compatible
BrandMax Triple Enzyme Cleaner	BrandMax	Compatible
Cetylite® Power Cleaner	Cetylite Industries Inc	Compatible
Crosstex® General Purpose Ultrasonic Cleaner	Crosstex	Compatible
Instrument Forte	Helvemed	Compatible
Maxizyme Tabs	Henry Schein	Compatible
Enzymax® Liquid	Hu-Friedy	Compatible
Enzymax® Powder	Hu-Friedy	Compatible
Enzymax Earth® Tablets	Hu-Friedy	Compatible
Patterson General Purpose Type1	Patterson	Compatible
Pearson Quality Super	Pearson	Compatible
Brite Sheild	Pemier	Compatible
ReSURGE™ Instrument Cleaning Solution	Sultan Healthcare	Compatible
Clean & Simple Tablets	Tuttnauer	Compatible
GZYME3	Germiphene Corp.	Compatible
Vigilance	Harry J. Bosworth Co.	Compatible
Citrizyme	Pascal International, Inc.	Compatible
Good Vibrations	Certol International LLC	Compatible
Resurge	Sultan Healthcare	Compatible
Empower Dual Enzymatic Detergent	Metrex Corporation	Compatible
Eco 4 Neutral All-Purpose Cleaner	Sonix4	Compatible
Z <sup>3</sup> ™ Multi-Purpose Enzymatic Detergent	Benco	Incompatible
neodisher® FA	Chemische Fabrik Dr. Weigert GmbH & Co. KG	Incompatible
neodisher® MA	Chemische Fabrik Dr. Weigert GmbH & Co. KG	Incompatible
neodisher® Mielclear	Chemische Fabrik Dr. Weigert GmbH & Co. KG	Incompatible
neodisher® N	Chemische Fabrik Dr. Weigert GmbH & Co. KG	Incompatible
Denta Zyme®	Ecolab	Incompatible
COEZYME™ Enzymatic Detergent	GC America	Incompatible
GZYMXE	Germiphene Corp	Incompatible
VibraKleen® E Tablets	Heraeis Kulzer	Incompatible
EZ-Zyme® All Purpose Cleaner	Integra Miltex	Incompatible
Multi-Purpose Enzyme	Patterson	Incompatible
Ultra Dose General Purpose Cleaner	L&R Manufacturing Company	Incompatible
Microbex Microzyme	Germiphene Corp.	Incompatible



AMERICAN EAGLE INSTRUMENTS® INC

better DENTISTRY by DESIGN®

XP® Technology

The Sharpen-Free Instruments™