


# Applicator 2.0

## Reinventing the Applicator

The Applicator is a completely new approach to the application of dental liquids and gels. Since it's not a brush and not made from fibers and glue, nothing contaminates the restoration or bonding materials. The Applicator picks up, holds, and releases the same exact dose every time, giving you precise control. With the Applicator, you finally have a modern, absolutely accurate tool that lets you work with precision and without contamination.

### The Applicator at a Glance:

- Non-shedding, will not contaminate the restoration site
- An exact amount of material dispensed, every time
- The first re-engineered applicator in 20 years



1:1

### A re-engineered solution

As an innovator and high-end manufacturer of dental tools, Transcendent has completely re-engineered a solution for the application of dental liquids and gels.

### No contamination

Transcendent recognized that today's complex bonding agents and dental materials should be applied with tools that don't contaminate bonding agents or the restoration site. For the last 20 years, all micro-applicators have been made the same way – fibers and glue. Fibers shed and contaminate. The Applicator combines a rigid handle with a soft pliable head. The head contains capillary canals that pick up and release precise doses of dental materials. The handle bends and holds its shape for tough-to-reach applications.

### More precision

Since the Applicator is precision molded, the capillary canals always hold and release the same amount of fluid every time you use one. Instead of painting with a brush, you can apply just the amount of liquid needed for the procedure, reducing waste and the risk of excess material in the restoration site. With a pliable head and no glue sphere, the Applicator flexes into tight spots and corners.

### To the highest standards

We make the Applicator with a manufacturing process that is totally automated and controlled to the highest standards. It is offered in color-coded soft and regular versions and in two sizes.

Get in touch with your dealer, or find out more at [www.appliquator.com](http://www.appliquator.com)



You do precision work.  
We create your tools.

transcendent®



Intro Kit REF 440





1:1

Four colours – two softnesses, two sizes –  
for flexible and precise application.

#### Regular

REF	QTY	Flexibility
441	400	normal 
442	400	soft 

#### Fine

REF	QTY	Flexibility
443	400	normal 
444	400	soft 

#### Intro Kit

REF	QTY	Kit
440	400	100 pieces each



#### Transcodent GmbH & Co. KG

Zum Kesselort 53  
24149 Kiel, Germany

Phone: +49 (0) 431 / 979 979 0  
Fax: +49 (0) 431 / 979 979 59

[info@transcodent.com](mailto:info@transcodent.com)  
[www.transcodent.com](http://www.transcodent.com)

You do precision work.  
We create your tools.

transcodent®