

Available Patterns

The Boge 513™ is available exclusively in XP Technology, in both fixed instruments and Quik-Tip® cone socket instruments.



AESB513XPX



AESB513XPZ



AMERICAN EAGLE INSTRUMENTS® INC
better DENTISTRY by DESIGN®

Presents



Boge 513™



AMERICAN EAGLE INSTRUMENTS® INC
better DENTISTRY by DESIGN®

www.am-eagle.com/Boge513 | 800.551.5172

Instrument Information

The Boge 513 is a modified anterior sickle/hoe scaler used to remove calculus in the dental prophylaxis process. This instrument is double ended, and is only available in sharpen-free XP Technology from American Eagle Instruments.

XP[®] Technology
The **Sharpen-Free** InstrumentsSM

Due to the angulation of this instrument, it would be difficult to sharpen, so producing it in XP is significantly beneficial to the clinician.



One end of the Boge 513 has an elongated shank and straight, slightly lipped cutting edge. The hoe's width was specifically designed to maintain stabilization during vertical and horizontal working strokes.



The alternate end of the instrument is a Nebraska 128 sickle scaler.

Clinical Application



The anterior hoe is designed to scale:

- marginal recession areas with exposed root surfaces
- both the upper and lower anterior intraoral regions
- areas where patients have fixed orthodontic wires
- posterior areas with root surface exposure



The Nebraska 128 sickle scaler is designed to scale:

- interproximal contact areas
- areas with crowding and/or recession



“The curvature of the blade is brilliant for lingual wires, and I have always loved my Nebraska 128 for challenging anteriors.”

- Ryan Rutar, RDH, MAM

About Emily Boge

Prior to beginning her career in administration, Emily practiced as a full-time clinical dental hygienist in Manchester, Iowa (2003-2015). She is the owner of Think Big Dental, an academic writing and speaking firm, and is the current President of the Iowa Dental Hygienists' Association (2016-2017). As an active member of her professional associations ADHA and AADH, Emily has advocated for the dental hygiene profession on many levels. Emily was presented with the 2015 Master Clinician Award from ADHA and was Iowa's IDHA 2012 Hygienist of the Year. Emily enjoys cooking, motorcycling, and spending time with her husband and two boys on their Century Farm north of Farley, Iowa.

