


ONE
FOR
ALL



GUIDED
BIOFILM
THERAPY®




BIOFILM MANAGEMENT WITH AIR-FLOW® PLUS POWDER




ABOVE
CEMENTO-
ENAMEL
JUNCTION

1. AIR-FLOW® SUPRAGINGIVAL



CEMENTO-
ENAMEL
JUNCTION
< 4mm

2. AIR-FLOW® SUB- AND SUPRAGINGIVAL



CEMENTO-
ENAMEL
JUNCTION
> 4mm

3. AIR-FLOW® SUBGINGIVAL

“I FEEL
GOOD”

EMS-DENTAL.COM

EMS
Electro Medical Systems SA
Chemin de la Vuarpillière 31
CH-1260 Nyon
Phone +41 22 99 44 700
Fax +41 22 99 44 701
contact@ems-ch.com

© EMS FA-569/EN - rev. A-01

ALL-IN-ONE
N°1
PROPHYLAXIS

ORIGINAL
METHOD

AIR-FLOW® POWDERS

GUIDED
BIOFILM
THERAPY®

EMS+

AIR-FLOW® POWDERS

ORIGINAL THERAPY AIR-FLOW®

→ SUB- AND SUPRAGINGIVAL



AIR-FLOW® POWDER PLUS


- SUB- AND SUPRAGINGIVAL - IN ONE STEP
- NO IRRITATION OF THE GUMS - PATIENTS FEEL GOOD AND FREE OF PAIN
- NO MEASURABLE IMPAIRMENT OF THE TOOTH SURFACE, EVEN IN RECALL PROPHYLAXIS
- Erythritol-based
- Average grain size ~14 µm
- 120 g bottle

CLINICAL APPLICATIONS

- REMOVAL OF BIOFILM AND MEDIUM/LIGHT DISCOLORATIONS WITH **OPTIMIZED PATIENT COMFORT**
- ELIMINATES SUPERFLUOUS POLISHING WITH POLISHING PASTE - THAT SAVES TIME

- REMOVAL OF BIOFILM FROM POCKETS **DOWN TO 4 MM** ON NATURAL TEETH AND IMPLANTS
- REMOVAL OF BIOFILM FROM DEEP PERIODONTAL POCKETS, ON NATURAL TEETH AND IMPLANTS
- WITH AIR-FLOW® PERIO HANDPIECE AND SOFT DISPOSABLE NOZZLE WITH DEPTH MARKS

AIR-FLOW® HANDPIECES



→ SUPRAGINGIVAL



AIR-FLOW® POWDER CLASSIC NEW FORMULA

- Sodium bicarbonate
- Average grain size ~40 µm
- 300 g bottle
- Four flavors: lemon, mint, cherry and neutral

- REMOVAL OF BIOFILM AND HEAVY DISCOLORATIONS WITH OPTIMIZED **PATIENT COMFORT** ABOVE CEMENTO-ENAMEL JUNCTION

