

NaviTip® FX® Tips

- One-of-a-kind brush cleans, scrubs, and irrigates simultaneously
- Rigid 30 ga cannula


Listed as an "EXCELLENT" product by a prominent independent research institute.¹

Designed for: Ultradent Citric Acid 20% Solution.

¹ Clinical Research Associates Newsletter, Volume 29, Issue 1, January 2005.

1452 - 17 mm 20pk

1454 - 25 mm 20pk



NaviTip® 31 ga with Sideport Irrigator Tips

- Double sideports deliver irrigants safely, minimizing the possibility of chemicals being expressed past the apex
- The world's smallest cannula navigates the most intricate canal spaces

Designed for: ChlorCid, Ultradent EDTA 18% Solution, Citric Acid 20% Solution, and Consepis.

31 ga (20pk)
5121 - 21 mm
5123 - 27 mm

31 ga (50pk)
5122 - 21 mm
5124 - 27 mm



NaviTip® Tips

- Provide controlled delivery to the apex
- Flexible, stainless steel cannulae easily navigate curved canals
- Available in 17, 21, 25, and 27 mm lengths and 29 or 30 ga cannulae

Designed for: EndoREZ, Ultradent EDTA 18% Solution, ChlorCid, ChlorCid V, Consepis, Consepis V, Ultradent Citric Acid 20% Solution, File-Eze, and UltraCal XS.

	29 ga (20pk)	30 ga (20pk)	29 ga (50pk)	30 ga (50pk)
17 mm, White	5112	1249	1378	1421
21 mm, Yellow	5113	1349	1374	1422
25 mm, Blue	5114	1250	1376	1423
27 mm, Green	5115	1354	1377	1424
Assorted	5116	1351	1379	3319



Quantity discounts available at 3 and 5 packs.

800.552.5512 ULTRADENT.COM

© 2016 Ultradent Products, Inc. All Rights Reserved.
68490.7 021516



Intraoral Tip

- Allows precise placement
- Attaches to dual-barrel mixing tips

Designed for: PermaFlo DC.

5922 - 20pk



Ultradent® Mixing Tip

- Mixes and delivers in one action

Designed for: UltraTemp REZ, UltraTemp, EndoREZ, and PermaFlo DC.

5920 - 20pk



Endo-Eze® Tips

- Work for a variety of endodontic procedures
- Come with flexible, strong cannulae in a variety of sizes
- Bend easily


Designed for: luting materials and air/water delivery.

Use with: TriAway Adapter, PermaFlo DC, and other Ultradent syringes.

20pk
345
346
347
348
350

100pk
1428
1429
1430
1431


Gauge - Diameter
18 ga - 0.0490";1.25 mm
19 ga - 0.0415";1.06 mm
20 ga - 0.0350";0.90 mm
22 ga - 0.0275";0.70 mm
30 ga - 0.0100";0.25 mm



ULTRADENT TIPS



DESIGNED TO DELIVER



ULTRADENT PRODUCTS, INC.

Improving Oral Health Globally



ULTRADENT PRODUCTS, INC.

800.552.5512 ULTRADENT.COM

The chemistries you use are different. Some are chemically activated, needing to be mixed immediately before delivery. Others have varying viscosities. Some work in pits and fissures, some inside canals, and some on enamel surfaces. Each chemistry you use is designed for a specific purpose. Shouldn't the same be true for your tips?

Ultradent makes tips **designed to deliver** each chemistry we create. Whether you're delivering a flowable composite, a viscous gel, or thick impression material, we make the perfect tip for the job. And since our tips are engineered on-site, we test each design to ensure it works perfectly with the chemistry it's intended for.



Our NaviTip® tips are treated to create flexible tips with a smooth, round end, ensuring they will easily snake through canals to deliver irrigants precisely where they are needed. This treatment can create a golden hue you may see on some of our NaviTips.

Check out our tips with **LOK-TITE®**

Luer Lock tips with Lok-Tite feature double threads that lock the tip into place for increased security, and wings for easy attachment and removal.



Black Micro® FX® Tip

- Accommodates various viscosities
- Flocked tip fans out to spread materials in a thin, uniform layer
- 22 ga tip

Designed for: PrimaDry, and PermaSeal.

1357 - 100pk 1434 - 500pk



Chromaclone® PVS Mixing Tip

- Enables direct delivery of impression materials
- Automixing, disposable, and color coded

Designed for: Chromaclone PVS.

2902 - Yellow 50pk 2903 - Pink 50pk 2904 - Green 50pk



Black Mini® Tip

- Dispenses large volumes
- Opaque plastic preserves the flow of light-cured materials

Designed for: Ultra-Blend plus, Ultradent LC Block-Out Resin, PermaFlo, TriAway Adapter, UltraTemp, Opalescence Boost, Ultradent Diamond Polish Mint, OpalDam, Opalescence Endo, OraSeal Caulking, and OraSeal Putty.

514 - 100pk 1433 - 500pk



Intraoral Impression Tip

- Allows precise placement
- Attaches to PVS Mixing tip

Designed for: Chromaclone PVS.

2906 - 50pk



Black Mini® Brush Tip

- Precise, controlled delivery of aqueous materials
- Tight, adjustable brush fibers minimize bubbles
- Unique to Ultradent

Designed for: Consepsis, Peak SE Primer, PermaFlo DC Primers, PermaQuick Primers, Seek/Sable Seek, Ultradent Silane, and Ultradent Universal Dentin Sealant.

193 - 20pk 1169 - 100pk 1432 - 500pk



SoftEZ® Tip

- Tip fibers provide visible, controlled delivery
- Brush fibers facilitate smooth application

Designed for: Enamelast.

4712 - 50pk 4711 - 100pk



Black Micro® Tip
Blue Micro® Tip

- Pinpoint precision
- Narrow cannulae accurately deliver materials

Black Micro designed for: Ultra-Blend plus, EtchArrest, OpalDam, UltraSeal XT plus, and UltraSeal XT hydro.
Blue Micro designed for: Ultra-Etch, Consepsis, Ultracare, and PrimaDry.

1085 - Black 100pk 1435 - Black 500pk
127 - Blue 100pk 1436 - Blue 500pk




White Mac™ Tip
White Mini™ Laser Tip

- Dispenses large volumes
- All-plastic delivery tips

Designed for: OraSeal Caulking, OraSeal Putty, Consepsis Scrub, Opalustre.

1361 - 100pk 1480 - 500pk 1247 - 20pk
White Mac White Mac White Mini Laser



Inspirational® Brush Tip

- Delivers viscous or filled materials smoothly via an internal helical channel and ridge
- Tight, adjustable brush fibers minimize bubbles

Designed for: Composite Wetting Resin, Peak Universal Bond, PermaQuick, Ultradent Porcelain Etch, Ultra-Etch, UltraSeal XT plus, and UltraSeal XT hydro.

123 - 100pk 1033-I - 500pk



Micro 20 ga FX® Tip

- Standard flowable delivery tip
- Flocked fibers fan out to evenly spread material over the entire surface

Designed for: Opalescence Boost.

1405 - 100pk 1406 - 500pk



Micro 20 ga Tip

- Large-gauge cannula enables consistent flow
- Standard flowable composite delivery tip

Designed for: PermaFlo, PermaFlo Purple, and OpalDam.

1252 - 100pk 1437 - 500pk



NanoTip® Tip

- Ideal for micro dentistry
- Small, strong, and rigid 30 ga cannula

Designed for: UltraSeal XT plus, Ultra-Etch, PermaFlo, microscopic restorative dentistry, and delivery of Vit-I-escence Colors for esthetic characterizations.

1240 - 20pk 1426 - 100pk



Metal Dento-Infusor® Tip

- Places hemostatic agents precisely, and effectively removes superficial coagulum
- Blunt, bent cannula with padded brush enables gentle pressure on the sulcus
- Ultradent's first tip, the "MDI" remains paramount for successful tissue management

Designed for: Astringedent, Astringedent X, ViscoStat, and ViscoStat Clear.

2559 - 100pk 2560 - 500pk




SST™ (Surgical Suction Tip)

- Ideal for delicate surgeries
- Large-diameter tip opening

Designed for: Luer Vacuum Adapter for small periodontic or endodontic surgeries and controlled suction of Opalescence Boost.

1248 - 20pk




Blue Mini® Dento-Infusor® Tip

- Offers the same tissue management benefits as the Metal Dento-Infusor tip
- Allows controlled flow of drop-sized quantities
- 1.20 mm tip

Designed for: Astringedent, Astringedent X, ViscoStat, and ViscoStat Clear.

1086 - 100pk 1440 - 500pk



Endo-Eze® Irrigator Tip

- Provides ideal reach without risk of expressing chemicals past the apex
- Comes with a flexible, blunt cannula with a unique, anti-obturating end
- 1.0" nonsterile and 1.25" sterile, 27 ga tip

Designed for: ChlorCid, ChlorCid V, and Ultradent 5 ml syringe.

207 - 20pk 202 - Sterile 20pk



Micro Capillary™ Tips

- Bright colors are easily identified against soft tissues
- The world's smallest molded tip (inside diameter: 0.008")

Designed for: Ultracare, perio materials, endodontics, and Luer Vacuum Adapter.

1120 - 5 mm 20pk 1121 - 10 mm 20pk 1802 - 5 mm 100pk



Capillary Tips

- Navigate canals and substantially minimize use of paper points
- Narrow, flexible taper accesses curved canals
- 0.014" opening, 0.36 mm tip
- 0.019" opening, 0.48 mm tip

Designed for: Luer Vacuum Adapter for moisture removal from endodontic canals.

341 - 0.014" 20pk 3099 - 0.014" 50pk
186 - 0.019" 20pk 1425 - 0.019" 50pk