Flexible Interproximal finishers with a



FitStrips are a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts, **Garrison Dental Solutions.**

One handy handle! FitStrips' simple, detachable handle makes a world of difference for both clinician ease of use and patient comfort. Your hand stays outside the mouth, granting a clear field of view and easier operation.

Do the Twist! Twist the color-coded barrel to adjust FitStrips' curvature and it's automatically locked in place. It doesn't get any easier than that!

FitStrips make your procedures easier!

- Composite restorations: contouring/ finishing all interproximal surfaces
- Orthodontics: IPR/tooth slenderizing
- Crown and bridge: cement removal and clean-up





Fit Strip Flexible I nterproximal finishers with a Twist



each grit (4 pcs)

Variety is the key to happiness!

One FitStrips starter kit contains everything you need:

- Single-sided diamond strips
- Double-sided diamond strips
- Multiple diamond grits
- Even serrated FitStrips are included!

ADDITIONAL KITS

All refills come in a package of 10 FitStrips with 2 handles.

Single Sided Kit 2 ea of all 4 grits plus 2 serrated	FPSK02
Double Side Kit 2 ea of all 4 grits plus 2 serrated	FPSK03

REFILLS

All refills come in a package of 10 FitStrips with 2 handles.

		Single Sided	Double Sided
Super fine 15 Micron 0.11 Thick	Ī	FPSXFSS	FPSXFDS
Fine 30 Micron 0.15 Thick		FPSFSS	FPSFDS
Medium 46 Micron 0.21 Thick		FPSMSS	FPSMDS
Coarse 90 Micron 0.30 Thick		FPSCSS	FPSCDS
Serrated 0.05 Thick		FPSSAW	

For more information go to www.garrisondental.com. Be sure to check out the website to see the full line of innovative Garrison products and our 6-month money-back guarantee.



To order: Call toll free 888.437.0032 or online at www.garrisondental.com

