



Style Italiano is a group of passionate clinicians and researchers, who are in constant development of materials and techniques for aesthetic and conservative dentistry. Their philosophy is based on developing solutions to make excellence in aestethic and functional dentistry attainable by all.

### **StyleItaliano**

founded by Dr. Walter Devoto and Prof. Angelo Putignano

For Polydentia, the collaboration with Style Italiano is more than just a source of inspiration. It is an opportunity to constantly renew its commitment to put excellence in Swiss manufacturing at the service of dentistry worldwide.

# Unica anterior STYLEITALIANÓ



Visit "Unica Community" to share Uniça anterior clinical cases

### Follow us #thinkunica



Polydentia SA via Cantonale 47 6805 Mezzovico-Vira Switzerland

phone: +41 91 946 29 48 fax: +41 91 946 32 03

info@polydentia.ch www.polydentia.ch



## the one solution for anterior restorations

## matrix system





#### ase of positioning a

3

2)

## **UNICa** anterior

powered by: STYLEITALIANO



"In today's conservative dentistry, the choice of products and the techniques we apply everyday must be predictable and costeffective, incorporating clear-cut advantages to our daily procedures. Unica is the solution to everyday anterior restorations, whether it be class III, IV, V, shape modifications, and/ or composite veneers in direct stratification technique. One effective matrix for everyday anterior aesthetic and functional results."

Dr. Jordi Manauta, DDS Head of StyleItaliano Idea Factory Sestri Levante, Italy



Unica is the simple and ideal matrix for anterior restorations such as class III, IV, V, direct stratification composite veneers, and shape modifications.

Thanks to its countoured shape, Unica adapts correctly to the different morphologies of anterior teeth and makes it possible to restore interproximal and cervical margins at once, even in presence of rubber-dam or gingival retraction cords, thus reducing chair time significantly. The placement wings allow fast and efficient matrix positioning.

Furthermore, Unica anterior matrix once positioned allows to easily visualize the final shape of the restoration, even before starting the procedure.



#### Class III restoration

In case of a Class III cavity where the proximal margin is missing, Unica is very useful as it allows the clinician to visualise the final shape even before proceeding with the restoration.

Dr. Giuseppe Chiodera Brescia, Italy









Composite veneers



Thanks to the new matrix Unica we have the possibility to restore at the same time the cervical and interproximal margins with significant advantages in terms of shape management and predictability of the restoration.

myWedge

Dr. Monaldo Saracinelli & Dr. Antonino Nicolò Italy







#### Class IV restoration



For this Class IV restoration, Unica helps to restore the cervical and interproximal profiles anatomically and in a very fast way.

Dr. Dan Lazar Oradea, Romania

polydentia





#### Shape modification

The unique design of Unica makes this matrix an essential and very efficient tool since it enables us to work both in the interproximal and cervical areas at the same time.

Dr. Carlos Fernández Villares – Madrid, Spain



