

# ZirClean™

## Restoration Cleaner

ZirClean is a cleaning gel designed for the non-abrasive cleaning of the bonding surfaces of zirconia (and other prosthetic restorations) after intraoral try-in. ZirClean helps to achieve reliable adhesive cementation results by removing the phosphate contamination of zirconia (as well as ceramic and metal restoration surfaces) that occur during try-in.\*

### Product Highlights


- ▶ Delivery system offers improved ease-of-use
- ▶ Easy placement and clean up
- ▶ Helps achieve reliable adhesive cementation results\*
- ▶ Contaminants such as saliva can be removed to help achieve better bond strength after try-in\*

### Clinical Significance

The use of ZirClean will help to achieve reliable cementation results with Zirconia restorations after intraoral try-in.\*



### Shear Bond Strength to Zirconia<sup>†</sup> (MPa) After Cleaning With...



<sup>†</sup>Bonding after saliva contamination using Ultrabond shear bond test using Z-Prime™ Plus and Duo-Link Universal™

### RELATED MATERIALS



ALL-BOND UNIVERSAL®



Z-Prime™ Plus

### ORDER INFORMATION

1 Syringe (5g), 2 Disposable Syringe Tip Pack (15pk),  
Instructions .....**B-7351P**

**800.247.3368 • 847.534.6000**

\*As compared to untreated samples. Data on file.  
Rx Only